






Reynard Premier Detergent & Disinfectant Wipes (RHS200, RHS205, RHS210, RHS216, RHS218), and **Neutral Detergent Wipes** (RHS201, RHS203, RHS206)
Compatibility Approval Confirmation by the following equipment and Device Manufacturers:

| | | Reynard Wipe Approved for Use | |
|---|---|--|--|
| | |  Premier Disinfectant Wipe Range |  Neutral Detergent Wipe Range |
| Company / Brand | Equipment or Device Type | | |
| Arjo Medical Beds | Enterprise 5000X, Enterprise 8000X, Citadel Patient Care Systems, and Citadel Plus | ✓ | ✓ |
| Baxter Healthcare | Evo IQ Infusion Pumps | ✓ | ✓ |
| BBraun | SPACE INFUSION PUMPS & ACCESSORIES • Syringe Perfusor® Space • Volumetric Infusomat® Space • Battery-Pack SP with WiFi • B. Braun SpaceStation • B. Braun SpaceStation MRI • Space GlucoseControl System (SGC System) • SpaceControl (SGC) | ✓ | ✓ |
| BD (Alaris) | VOLUMETRIC SYRINGE PUMPS • BD Alaris™ neXus GP Volumetric Pump | ✓ | ✓ |
| Canon Medical Systems | ULTRASOUND SYSTEMS & TRANSDUCERS (excepting: Canon Transesophageal and Laparoscopic probes as these specialised probes should go through a sterilization process) | ✓ | ✓ |
| Crescent Healthcare | PATIENT TROLLEY MATTRESSES • Premiflex • Standard Combination | ✓ | ✓ |
| Cubro | HOSPITAL BEDS & ACCESSORIES • Alrick bariatric • Bock • Cubro • Liberty Healthcare • Wissner-bosserhoff | ✓ | ✓ |
| Daniels Health | REUSABLE SHARPS CONTAINERS • SharpSmart Sharps Containers • Trolleys & Accessories | ✓ | ✓ |
| Dr Mach, Hospital Lighting Systems | Full range of hospital grade and medical lights. | ✓ | ✓ |
| | | NB: Approved for these models only: LED series: 2,3,5, 6 and 8. | |

Reynard Premier Detergent & Disinfectant Wipes (RHS200, RHS205, RHS210, RHS216, RHS218), and **Neutral Detergent Wipes** (RHS201, RHS203, RHS206)
Compatibility Approval Confirmation by the following equipment and Device Manufacturers:

| | | Reynard Wipe Approved for Use | |
|---------------------------|---|---|---|
| Company / Brand | Equipment or Device Type |  |  |
| | | Premier Disinfectant Wipe Range | Neutral Detergent Wipe Range |
| Dräger | PATIENT MONITORING SYSTEMS • M540 ANAESTHESIA MACHINES • Zeus IE | | ✓ |
| Fresenius Kabi | AGILIA INFUSION PUMPS & ACCESSORIES • VP MC • SP MC • SP TIVA • SP PCA • Link • Duo | ✓ | ✓ |
| FUJIFILM SonoSite | ULTRASOUND SYSTEMS & TRANSDUCERS -please email sales@reynardhealth.com to obtain the full list of compatible Systems and transducers) | ✓ | ✓ |
| GE Healthcare | ULTRASOUND SYSTEMS & TRANSDUCERS -please email sales@reynardhealth.com to obtain the full list of compatible Systems and transducers) | ✓ | ✓ |
| Hillrom | HOSPITAL BEDS AND ACCESSORIES | ✓ | ✓ |
| Howard Wright | HOSPITAL BEDS AND MATTRESSES | ✓ | ✓ |
| Mindray | PATIENT MONITORING SYSTEMS • BeneVision N – 12 • BeneVision N – 15 • BeneVision N – 17 | ✓ | ✓ |
| Pharmaco (NZ) Ltd | BLOOD SALVAGE • Fresenius Kabi Continuous Auto Transfusion System • Hologic Fetal Fibronectin Analyser | ✓ | ✓ |
| Philips | ULTRASOUND SYSTEMS & TRANSDUCERS -please email sales@reynardhealth.com to obtain the full list of compatible Systems and transducers) DIGITAL DIAGNOSTIC RADIOLOGY SYSTEMS | | ✓ |
| Siemens Healthcare | ANGIOGRAPHY SYSTEMS • ARTIS icono • ARTIS pheno systems ULTRASOUND SYSTEMS & TRANSDUCERS • P500 System • Sequoia System • Redwood System • Juniper System Also approved for use for the majority of Siemens Transducers | ✓ NB: Premier Wipes are not approved for use on the Monitors of the Sequoia, Redwood and Juniper Systems | ✓ |

Reynard Premier Detergent & Disinfectant Wipes (RHS200, RHS205, RHS210, RHS216, RHS218), and **Neutral Detergent Wipes** (RHS201, RHS203, RHS206)
Compatibility Approval Confirmation by the following equipment and Device Manufacturers:

| | | Reynard Wipe Approved for Use | |
|-----------------------|--|---|--|
| Company / Brand | Equipment or Device Type |  Premier Disinfectant Wipe Range |  Neutral Detergent Wipe Range |
| Smiths Medical | CADD® INFUSION PUMPS & ACCESSORIES • CADD®- Solis VIP Ambulatory - Infusion Pump • CADD®- Solis Ambulatory Infusion - Pump with Wireless Communications • CADD-Legacy® PCA Pump • CADD-Legacy® PLUS Pump • CADD® Pump Lockbox | | ✓ |
| Toshiba | ULTRASOUND SYSTEMS & TRANSDUCERS (excepting: Toshiba Transesophageal and Laparoscopic probes as these specialised probes should go through a sterilization process) | ✓ | ✓ |
| Welch Allyn | PATIENT MONITORING SYSTEMS • Connex Spot Monitor (CSM) | ✓ NB: Premier wipes are only compatible for the main housings only. The LCD displays and any accessory components (including Braun thermometer handles) should follow cleaning instructions as provided in the respective Welch Allyn Directions for Use | ✓ |

Reynard Health Supplies. Instrument Compatibility Approval Update 24/8/2021.
 Issue Number 8.

NB: If there are any other devices and/or medical equipment used within your healthcare setting that are not represented on this approval list. Please email Reynard directly at sales@reynardhealth.com with the equipment details, and we will arrange to have compatibility testing conducted with the relevant manufacturer. Thank you.



Reynard is reassurance.



Reynard
 Health Supplies

New Zealand
 PO Box 8470
 Havelock North
Sales@reynardhealth.com
 0508 494 737

Australia
 44 Hampden Road
 Artarmon, NSW, 2064
enquiries@reynardhealth.com.au
 1300 667 199

UK
 Boston House, Grove Business Park
 Downsview Road, Wantage
 OXFORDSHIRE, OX12 9FF
christine@reynardhealth.com
 +44 (0)845 263 7335

International
 PO Box 8470
 Havelock North
Sales@reynardhealth.com
 +64 (0)6 650 0709