

Reynard Microfibre Mop Heads are a fast & easy way to clean your floors.

Microfibre Mop Heads are designed to deliver a high quality clean and reduce the potential for cross-contamination. Microfibre technology is proven to eliminate 99% of all dirt, dust and pathogens. The unique fibres do all the cleaning thus reducing the mechanical action needed. Can be used with or without steam, chemicals and other cleansers.

Specifications


- Tested by consultant Microbiologist
- Removed 99.9% of bacteria

Points to remember


- Always start from cleanest area to dirty mopping in a S shape pattern
- Always dispose of mops if soiled into appropriate receptacle

80ml	400ml	800ml	1.2 Litres	1.6 Litres	2 Litres
1 Mop	5 Mops	10 Mops	15 Mops	20 Mops	25 Mops


Step 1: Place the required number of mop heads into a bucket



Step 2: Pre-soak the mop heads with the required amount of water as per the table above. Turn the mop heads over if required when you first add the water.



Step 3: Wait for the water to fully absorb before using. Drain off excess water.



Step 4: Start from cleanest to dirtiest area. Mop in a S pattern.



CODE	COLOUR	UOM
RHS318	White	25 mop heads per inner bag 4 inner bags per carton

Reynard is reassurance.



New Zealand
PO Box 8470
Havelock North
Sales@reynardhealth.com
0508 494 737

Australia
44 Hampden Road
Artarmon, NSW, 2064
enquiries@reynardhealth.com.au
1300 667 199

UK
White Horse Business Park,
Standford in the Vale SN7 8NY
www.reynardhealth.co.uk
uk@reynardhealth.com
+44 (0)845 263 7335

International
PO Box 8470
Havelock North
Sales@reynardhealth.com
+64 (0)6 650 0709