

Reynard Microfibre Mop Heads.

Chemical free cleaning.

Designed to deliver a high quality clean while reducing the potential for cross contamination. Proven to eliminate 99.9% of all dirt, dust, and pathogens without the use of harsh chemicals. Ideal for use wet, dry, or used with steam for the most effective clean.

Points to remember:

- Mop from clean to dirty
- Always mop in an 'S' shape pattern
- Dispose of mops if soiled into appropriate receptacle.




80ml	400ml	800ml	1.2 Litres	1.6 Litres	2 Litres
1 Mop	5 Mops	10 Mops	15 Mops	20 Mops	25 Mops

Step 1: Place the required number of mop heads into a bucket.

Step 2: Pre-soak the mop heads with the required amount of water as per the table above. Turn the mop heads over if required when you first add the water.

Step 3: Wait for the water to fully absorb before using. Drain off excess water.

Step 4: Start from the cleanest to the dirtiest area. Mop in an 'S' pattern.



Product Code	RHS318
Heads per bag	25
Bags per carton	4



Reynard is reassurance.

